

PROTOCOL FOR *Ukko*;Hk2^{+/-} MOUSE GENOTYPING

Procedure

Genotyping of offspring from Ukko;Hk2^{+/-} breeding colony is based on PCR. Run all three primers in the same reaction.

PCR primers

5' forward primer (Hk2KO-F2) - 5' **act ctc ctg ccg ccc tgc** 3'

3' reverse primer (Hk2KO-R1) - 5' **ccc ctc atc gcc acc gc** 3'

3' reverse primer (Neo-R1) - 5' **gtg ccc agt cat agc cga ata gc** 3'

PCR profile; programme UKKO

95 °C, 6 min

94 °C, 30 s

40 cycles

56 °C, 30 s

72 °C, 1 min

72 °C, 10 min

 $4^{\circ}\text{C}, \infty$

PCR mix

10 x PCR Gold buffer (Perkin Elmer)	3.0 μ l
MgCl ₂ (25 mM)	1.5 μ l
dNTPs (10 mM)	0.5 μ l
Hk2KO-F2 (20 μ M)	1.0 μ l
Hk2KO-R1 (20 μ M)	0.5 μ l
Neo-R1 (20 μ M)	0.5 μ l
AmpliTaq Gold (5 U/ μ l)	0.2 μ l
DNA template (~ 0.5 μ g tail DNA)	2.0 μ l
H ₂ O	<u>20.8 μl</u>
	30 μ l

Post-PCR analysis

Load 12 μ l of the PCR reaction on a 1% agarose gel.

Expected pattern; for wt 239 bp and for mutant 432 bp.