


GENOTYPING ROSA-DICRE MICE

I. GENOTYPING ROSA-DICRE MICE USING PRIMERS SPECIFIC FOR THE DICRE CONSTRUCT (PLOS ONE, 2007:2(12), E1355)

PRINCIPLE: Amplification of the Cre59-SV40 polyA junction



PCR PROCEDURE

Primers :

| Target | Primer pair | Amplicon size (bp) |
|--------------------------------|---|--------------------|
| Rosa-DiCre DNA (iCre59-SV40pA) | Pr236 <i>acctctgatgaagtcaggaagaa</i> Pr264 <i>ccacaactagaatgcagtga</i> | 411 bp |

Mix:

| | Volumes (µl) for 1 reaction | Volumes (µl) for 25 reactions |
|--|-----------------------------|-------------------------------|
| 10x Taq Platinum buffer (w/o MgCl ₂) | 2 | 50 |
| dNTP, 10 mM stock | 0,4 | 10 |
| Primers, 10 pmol/µl stock | (0,4 + 0,4) | (10 +10) |
| MgCl ₂ 50 mM | 1,2 (3mM final) | 30 |
| BSA, 10 mg/ml (Biolabs) | 1 | 25 |
| Taq Platinum (Invitrogen, ref 10966034) | 0,08 | 2 |
| H ₂ O, qsp 19 µl | 13,52 | 340 |


Put together at RT, distribute 19 µl in each PCR tube.

Add 1 µl target DNA prepared by the HotShot procedure (Truett et al., Biotechniques 29:52-54, 2000).

Do PCR, using 57°C for annealing, 30 secs elongation time and 35 cycles.

Run a 2% agarose gel loaded with 10 µl of the reaction.

EXPECTED RESULTS



Exemple with three different annealing temperatures, with DiCre DNA (positive band) and a negative control (Cre2ERT2 DNA)

POSSIBLE ALTERNATIVES

Use the same PCR mix and conditions, but using either of the following primer pairs:

| Target | Primer pairs | Amplicon size (bp) |
|-----------------------|---|--------------------|
| FRB-Cre60 junction #1 | Pr109 <i>taatggaggccaagagtg</i> Pr235 <i>gcatccacattctcctttctg</i> | 344 bp |
| FRB-Cre60 junction #2 | Pr280 <i>ttggggaaggaacgtgaaagg</i> Pr279 <i>aggtgctgttgatggtctca</i> | 355 bp |

II. GENOTYPING ROSA-DICRE MICE USING THE GENERIC PRIMERS OF P. SORIANO, DETECTING THAT THERE IS AN INSERT (ANY) WITHIN THE ROSA26 LOCUS.

NOTES:

1. This procedure allows the **identification of homozygous vs. heterozygous vs. WT animals**.
2. It works when insertion has been done using P.Soriano's targeting construct.
3. The conclusion is based on the absence of the 603 bp band corresponding to the WT locus => don't forget to have a WT or a heterozygous control (to test that PCR works OK)

Perform 2 separate PCR reactions (multiplexing is hazardous due to competition between primers), using the following primer pairs:

| Target | Primer pairs | Amplicon size (bp) |
|--------------------------|---|--------------------|
| Wild-type Rosa26 locus | Pr306 <i>ggagcgggagaaatggatag</i> Pr307 <i>aaagtcgctctgagttgttat</i> | 603 |
| Rosa26 locus with insert | Pr307 (as above) Pr305 <i>gccaagagtttgcctcaacc</i> | 325 |

Use PCR mix and conditions as above, but set annealing at 60°C.