

## Noggin genotyping by PCR Protocol adapted on 10/10/12

The Noggin genotyping is done by 2 PCRs, 1 PCR to amplify the wild type gene and 1 to amplify the knockout gene.

Prepare the mix on ice.

PCR mix per sample.

	WT ( $\mu$ l)	KO ( $\mu$ l)
H <sub>2</sub> O (Fresh Milli-Q)	6.5	6.5
Buffer	1	1
MgCl <sub>2</sub>	0.5	0.5
dNTP	1	1
primers 210 x 211	1	
primers 210 x 212	-	1
Enzyme (Taq-Pol.)	0.15	0.15
Template DNA	1,5	1,5
Betaine 5M	2	2

### Mix by pipetting

PCR program: 2-step PCR.

95 °C - 2 min	
95 °C - 20 seconds	} 32 cycles
72 °C - 1 min	
72 °C - 3 min	
4 °C - $\infty$	

The number of cycles may be increased depending on the quality of DNA (for instance not isolated by the machine)

Expected size:

Around 200 basepairs for the WT PCR (210 x 211)

Around 170 basepairs for the KO PCR (210 x 212)

Results:

Only the WT PCR amplified the fragment → mouse is a wildtype

Only the KO PCR amplified the fragment → mouse is a knockout

Both PCRs amplified the fragments → mouse is a heterozygote

Remarks:

Use buffer, MgCl<sub>2</sub> and enzyme from the same kit (Lot#).

Don't forget the control samples. There is a special plate with control samples for all PCRs in fridge 3 in a box (which normally contains 96-Well Optical Reaction Plates for purifying DNA).

PRIMERS:

- primers for WT pcr

F: gcatggagcgctgccccagc

R: gagcagcgagcgagcagcagcg

- primers for KO pcr

F: gcatggagcgctgccccagc (same F primer as for WT pcr)

R: aagggcgatcggtgcgggcc

PRODUCTS:

- Buffer, MgCl<sub>2</sub> and Taq polymerase: Bioline Biotaq polymerase 2500 units; cat nr BIO-21060

- dNTP: Life technologies 10mM dntp mix; cat nr 18427-013

- Betain: Sigma-Aldrich 5M betain; cat nr B0300-1VL