

Cre genotyping

The Cre genotyping is done by 1 PCR which amplifies the Cre gene. But since matings only take place between heterozygotes and wild types, progeny can only be heterozygotes or wild types.

| Name | Oligo Sequence |
|-------------------|--------------------------|
| PT178-Cre forward | TCCAATTTACTGACCGTACACCAA |
| PT179-Cre reverse | CCTGATCCTGGCAATTTCCGGCTA |

Prepare the mix on ice.

PCR mix per sample.

| | (μ l) |
|----------------------------------|------------|
| H ₂ O (Fresh Milli-Q) | 6.5 |
| Buffer | 1 |
| MgCl ₂ | 0.5 |
| dNTP | 1 |
| primers 178 x 179 | 1 |
| Enzyme (Taq-Pol.) | 0.1 |
| Template DNA | 1 |

Mix by pipetting

PCR program: 3-step PCR.

| | |
|--------------------|-------------|
| 95 °C - 1 minute | |
| 95 °C - 10 seconds | } 30 cycles |
| 60 °C - 10 sec | |
| 72 °C - 30 sec | |
| 72 °C - 3 min | |
| 4 °C - ∞ | |

Expected size:

550 basepairs

Results:

No fragment → mouse is a wildtype

Fragment of 550 basepairs → mouse is a heterozygote

Remarks:

Use buffer, MgCl₂ and enzyme from the same kit (Lot#).

Buffer, MgCl₂, dNTP and primers are stored in fridge 3 in Eric's PCR Equipment Box, the Taq-Polymerase is stored in freezer 5, top drawer on the left side in Eric's Box.

Don't forget the control samples. There is a special plate with control samples for all PCRs in fridge 3 in a box (which normally contains 96-Well Optical Reaction Plates for purifying DNA). This box is called **DNA Eric**.