



Genotyping protocol for $\Delta 9$ mice:



Targeted genomic locus of mouse *Cyld*



Primers used

FWD1: 5'-GATGGCTCTTGTCACCACTT-3'

F5: 5'- GCA GGC TGT ACA GAT GGA AC- 3'

F6: 5'-CGT GAA CAG ATG TGA TGA AGG-3'

R6: 5'-CTA CCA TCC CTG CTA ACC AC-3'

For the detection of WT allele, we use primer combination F6 x R6. For the detection of KO allele, we use either primer combination FWD1 x F6 x R6 (new protocol) or primer combination F5/ R6 (old protocol)

WT allele Protocol

F6: 5'-CGT GAA CAG ATG TGA TGA AGG-3'

R6: 5'-CTA CCA TCC CTG CTA ACC AC-3'

PCR program

- 1) T=94°C 3 min
- 2) T=94°C 45 sec
- 3) T=54°C 45 sec (depending on the PCR machine, range 52-56°C)
- 4) T=72°C 1min 30 sec
- 5) Go to step 2 35x
- 6) T=72°C 7 min
- 7) T=16°C forever

Expected size: 350 bp (WT allele), 450bp (Floxed allele). For $\Delta 9$ heterozygotes we detect only the 350bp band, whereas for homozygote $\Delta 9 / \Delta 9$ mice F6XR6 combination doesn't give any product.

Reaction


Template: we use 2 μ l DNA from 1/20 dilution of 30 μ l genomic DNA preparation from 2mm mouse ear punch. (Each mouse ear punch yields typically 2-3 μ g genomic DNA)

2,5 μ l dNTPs (2 mM stock)5 μ l 5x reaction buffer1,5 μ l MgCl₂ (stock 25 mM)1 μ l primer F6 (10 pmol/ μ l)1 μ l primer R6 (10 pmol/ μ l)0,25 μ l TaqH₂O up to 25 μ l

Analyze PCR product in 2% agarose gel.

Controls: C57BL/6 genomic DNA, Fl9/wt genomic DNA

Expected image for Floxed/WT bands, using F6-R6 primer combination.



KO allele genotyping protocol, using primer combination FWD1 x F6 x R6
(NEW protocol)

FWD1: 5'-GATGGCTCTTGTCACCACTT-3'

F6: 5'-CGT GAA CAG ATG TGA TGA AGG-3'

R6: 5'-CTA CCA TCC CTG CTA ACC AC-3'

PCR program (Biorad machine)

- 1) T=94°C 3 min
- 2) T=94°C 45 sec
- 3) T=57°C 45 sec
- 4) T=72°C 1min 30 sec
- 5) Go to step 2 35x
- 6) T=72°C 3 min
- 7) T=4°C forever

Expected size: 250 bp for KO allele and 350bp for WT

Reaction

Template: we use 2 μ l DNA from 1/20 dilution of 30 μ l genomic DNA preparation from 2mm mouse ear punch. (Each mouse ear punch yields typically 2-3 μ g genomic DNA)

2,5 μ l dNTPs (2 mM stock)

5 μ l 5x reaction buffer

1,5 μ l MgCl₂ (stock 25 mM)

1 μ l primer F6 (10 pmol/ μ l)

1 μ l primer FWD1 (10 pmol/ μ l)

2 μ l primer R6 (10 pmol/ μ l)

0,25 μ l Taq

1 μ l DMSO

H₂O up to 25 μ l


Analyze PCR product in 2% agarose gel.

Controls: C57BL/6 genomic DNA, KO MEF DNA ($\Delta 9/\Delta 9$) and or $\Delta 9/+$ genomic DNA.

There is a preference for KO allele using the triple primer combination, when having heterozygotes.

Expected image for WT/KO bands using F6-FWD1-R6 primer combination.

PCR for screening of $\Delta 9$ and WT alleles, using a set of 3 primers (F6-R6-KO1)



KO allele genotyping protocol, using primer combination F5 x R6
(OLD protocol)

F5: 5'- GCA GGC TGT ACA GAT GGA AC- 3'

R6: 5'-CTA CCA TCC CTG CTA ACC AC-3'

PCR program

1. T=94°C 3 min
2. T=94°C 45 sec
3. T=52°C 45 sec
4. T=72°C 90 sec
5. Go to 2 35x
6. T=72°C for 3 min
7. Hold at 16°C

Expected size: 700 bp

Reaction

Template: we use 1 μ l DNA from 50 μ l genomic DNA preparation from mouse tails. (Dissolve @ 56°C for 1-2h before use.)

- 2,5 μ l dNTPs (2 mM stock)
- 5 μ l 5x reaction buffer
- 1,5 μ l MgCl₂ (stock 25 mM)
- 1 μ l primer F5 (10 pmol/ μ l)
- 1 μ l primer R6 (10 pmol/ μ l)
- 0,25 μ l Taq
- H₂O up to 25 μ l

Analyze PCR product in 2% agarose gel.

Controls: C57BL/6 genomic DNA, KO MEF DNA ($\Delta 9/\Delta 9$) or $\Delta 9/+$ genomic DNA.

Expected image for WT/KO bands, using F5-R6 primer combination.

