

Dr Vassilis Aidinis laboratory  
Immunology Research Department  
Alexander Fleming Biomedical Sciences Research Center  
34 Fleming Street,  
Vari, 16672  
Greece  
Tel. +30 210 9656310 (ext. 159)  
Fax +30 210 9653934  
Email: aidinis@fleming.gr

## **Genotyping protocol for lysophosphatidic acid receptor 1 (LPAR1) (MGI ID: 108429) lacZ reporter mice.**

Protocol edited by Evangelos Dioletis on 22-10-2015.

Information presented here applies to the lacZ-tagged null allele ( $Lpar1^{tm1b(EUCOMM)Wtsi}$ ), which was derived from the  $Lpar1^{tm1a(EUCOMM)Wtsi}$  allele via Cre-mediated recombination .

### Genotyping primers

Fwd primer (727): 5'- CGGTCGCTACCATTACCACT -3' (20 bp,  $T_m=60.5^\circ\text{C}$ )  
Rev primer (728): 5'- ACTGATGGCGAGCTCAGACC -3' (20 bp,  $T_m=62.5^\circ\text{C}$ )

Primers 727 and 728 bind mutant gDNA only.  
Expected PCR product length: 380 bp

The current protocol does not discriminate between homo- and heterozygous mice.

Typical pipetting scheme for single reactions (in  $\mu\text{l}$ ):


ddH <sub>2</sub> O	12.4
10x buffer*	2
MgCl <sub>2</sub> (25mM)	1.2
dNTPs (2.5 mM each)	1
Fwd primer (5 $\mu\text{M}$ )	1 (727)
Rev primer 1 (5 $\mu\text{M}$ )	1 (728)
DNA template (50-100 ng/ $\mu\text{l}$ )	1
Taq DNA polymerase	0.4
<b>Total volume (<math>\mu\text{l}</math>)</b>	<b>20</b>

\*10x PCR Buffer: 500 mM KCl, 100mM Tris-HCl (pH 9.0 at 25°C), 1% Triton X-100.

Thermocycler conditions:

94°C for 5 min  
94°C for 30 s  
58°C for 30 s } x35  
72°C for 45 s  
72°C for 5 min

Example of genotyping LPAR1 lacZ reporter (*Lpar1*<sup>tm1b(EUCOMM)Wtsi</sup>) mice.



2% agarose gel electrophoresis in 1x TBE  
35 min at 120V (6.7V/cm).

Controls

M : Marker (*Pst*I-cut lambda DNA)  
H<sub>2</sub>O : water instead of gDNA template (-ve control)  
+ve : gDNA from LPAR1 tm1b mouse (+ve control)

## ***Pst*I-cut lambda DNA**




Figure taken from: [http://www.taq-dna.com/phage-lambda-dna-psti-digest-ready-to-use-\\_144.html](http://www.taq-dna.com/phage-lambda-dna-psti-digest-ready-to-use-_144.html)