

Genotyping protocol

General information:

Strain name	B6-Sdhd-KO
-------------	------------

Primers for PCR 1 and PCR 2:

Name	Sequence	Primer type
501	GGAAGGCTCCAAGGGTGCAG	gene-specific
502	CACATACACGCAGGCACTGG	gene-specific
505	TCTCTGCAGTGGCCAAGG	gene-specific
504	ATCTTGTTCAATGGCCGATCCC	selection cassette

In case more than two primers are introduced, please indicate how they should be combined:

	Forward primer	Reverse primer
Wild type	501	502
Null allele	505	504

Reaction mix for wt allele detection (PCR 1):

DNA	2	μ l
DreamTaq Green PCR Master Mix (2x)	12.5	μ l
501 primer 20 μ M	1	μ l
502 primer 20 μ M	1	μ l
H ₂ O	8.5	μ l
Final volume	25	μ l

PCR 1 program:

95	°C	3	min	X35
95	°C	30	sec	
65	°C	30	sec	
72	°C	60	sec	
72	°C	10	min	

Expected fragment size:

wt	370	bp
----	-----	----

Comments/Additional information:

Null allele detection requires a different PCR reaction (see below)

Reaction mix for null allele detection (PCR 2):

DNA	2	μl
DreamTaq Green PCR Master Mix (2x)	12.5	μl
505 primer 20 μM	1	μl
504 primer 20 μM	1	μl
H ₂ O	8.5	μl
Final volume	25	μl

PCR 2 program:

95	°C	3	min	X35
95	°C	30	sec	
65	°C	30	sec	
72	°C	60	sec	
72	°C	10	min	

Expected fragment size:

Null allele	430	bp

Comments/Additional information:

The targeting construct contains a neo cassette. If *Sdhd*^{+/-} mice are crossed with CRE-expressing lines, neo cassette can be excised from the null allele and PCR-3 will be needed for genotyping (see below)

Primers for PCR 3:

Name	Sequence	Primer type
505	TCTCTGCAGTGGCCAAGG	gene-specific
506	GAGTCTATTTCAAACATGTAAGCACC	gene-specific

Reaction mix for « excised » null allele detection (PCR 3):

DNA	2	μl
DreamTaq Green PCR Master Mix (2x)	12.5	μl
505 primer 20 μM	1	μl
506 primer 20 μM	1	μl
H ₂ O	8.5	μl
Final volume	25	μl

PCR 3 program:

95	°C	3	min	X35
95	°C	30	sec	
58	°C	30	sec	
72	°C	60	sec	
72	°C	10	min	

Expected fragment size:

« Excised »	300	bp
Null allele		