

Name of Mouse model or mutation:

CXCR4-CRE-ERT2-EM2-B6N

Description:

ERT2 cassette knock in to the existing *Cxcr4*-CRE-EM1-B6N allele made by CRISPR/Cas9 gene editing.

Type of mutation:

Knock-in: ERT2 cassette

Delivery method:

Pronuclear injection into 1-cell stage embryo.

Genetic Background:

C57BL/6NTac

Nuclease:

Cas9 mRNA

sgRNAs:

Protospacer sequence	PAM sequence
aggggcaatggtgcgctgc	TGG
cttaatcgccatctccagc	AGG
tcctatactattgaagaat	AGG

IssDNA donor sequence template (5'-3', NOTE donor delivered to embryo will be reverse complement (Codner et al., 2018: Figure 1)):

LOCUS *Cxcr4_Cre_ERT2_D* 1373 bp DNA linear 08-DEC-2021

FEATURES Location/Qualifiers

PCR_primer complement(3..22)
 /dnas_title="p_Cxcr4_2cuts_seqR4"
 /note=""
 /note=""
 /current=0

PCR_primer 89..108
 /note="Cre_F3"

```

        /current=0
misc_feature 1..200
        /note="200 nt 5'HA"
misc_feature 184..200
        /note="Cxcr4_CRE-ERT2_5_1"
misc_feature complement(200..217)
        /note="Cxcr4_CRE-ERT2_5_2 (OTs in B6N)"
misc_feature 197..199
        /note="SMs to kill the activity of guides 5_1 and 5_2:
        cgc>AgA R=R usage from 9.4 to 12.1"
misc_feature 1..217
        /note="CRE"
misc_feature 227..1171
        /note="ERT2"
misc_feature 218..226
        /note="Linking sequence as per Adipoq Cre-ERT2"
misc_feature 1175..1185
        /note="Remains of F3"
PCR_primer 1258..1280
        /dnas_title="p_Cxcr4_2cuts_seqF2"
        /note=""
        /note=""
        /current=0
PCR_primer complement(1326..1349)
        /dnas_title="p_Cxcr4_2cuts_seqR5"
        /note=""
        /note=""
        /current=0
misc_feature 1186..1373
        /note="Sv40 polyA site"
misc_feature 1175..1373
        /note="200 nt 3'HA"
misc_feature complement(1350..1373)
        /dnas_title="Cxcr4-ERT2_RNA_R"
        /note="Cxcr4-CRE-ERT2_D2_RNA_R"
source 1..1373
        /dnas_title="Cxcr4_Cre_ERT2_Donor2"

```

ORIGIN

```

1 gatttacggc gctaaggatg actctgggtca gagatacctg gcctgggtctg gacacagtg
61 ccgtgtcgga gccgcgcgag atatggcccg cgctggagtt tcaataccgg agatcatgca
121 agctgggtggc tggaccaatg taaatattgt catgaactat atccgtaacc tggatagtg

```

181 aacaggggca atgggAgAc tgctggaaga tggcgatctc gagccatctg ctggagacat
241 gagagctgcc aaccttggc caagcccgt catgatcaaa cgctctaaga agaacagcct
301 ggccttgcc ctgacggccg accagatggt cagtgccttg ttggatgctg agcccccat
361 actctattcc gagtatgac ctaccagacc cttcagtga gcttcgatga tgggcttact
421 gaccaacctg gcagacaggg agctgggtca catgatcaac tgggcgaaga ggggtgccagg
481 ctttgggat ttgaccctcc atgatcaggt ccaccttcta gaatgtgcct ggctagagat
541 cctgatgatt ggtctcgtct ggcgctccat ggagcacca gtgaagctac tgtttgctcc
601 taactgctc ttggacagga accagggaaa atgtgtagag ggcaggtgg agatcttca
661 catgctgctg gctacatcat ctcggttccg catgatgaat ctgcaggag aggagttgt
721 gtgcctcaaa tctattatt tgcttaattc tggagtgtac acatttctgt ccagcacct
781 gaagtctctg gaagagaagg accatatcca ccgagtcctg gacaagatca cagacactt
841 gatccacctg atggccaagg caggcctgac cctgcagcag cagcaccagc ggctggcca
901 gctcctctc atcctctccc acatcaggca catgagtaac aaaggcatgg agcatctgta
961 cagcatgaag tgcaagaacg tgggtcccct ctatgacctg ctgctggagg cggcggacgc
1021 ccaccgcta catgcgcca ctagccgtgg aggggcatcc gtggaggaga cggaccaaag
1081 cacttggcc actcgggct ctacttcac gattccttg caaaagtatt acatcacggg
1141 ggaggcagag ggtttcctg ccacagctg atgaATAGGA ACTTCgatca taatcagcca
1201 tatcacatct gtagaggttt tacttgcttt aaaaaactc ccacacctc cctgaacct
1261 gaaacataaa atgaatgcaa ttgttgtgt taactgttt attgcagctt ataatggtta
1321 caataaagc aatagcatca caatttcac aataaagca ttttttcac tgc

//

Pronuclear Microinjection mixes:

Microinjection buffer (MIB; 10 mM Tris–HCl, 0.1 mM EDTA, 100 mM NaCl, pH7.5) was prepared and filtered through a 2 nm filter and autoclaved. Cas9 mRNA, sgRNAs and ssODNs were diluted and mixed in MIB to the working concentrations of 100 ng/μl, 50 ng/μl each and 50 ng/μl, respectively. Injected embryos were re-implanted in CD1 pseudo-pregnant females. Host females were allowed to litter and rear F₀ progeny.

Sequence details

Cxcr4-CRE-EM1-B6N sequence (allele targeted)

ccacgaccaagtgacagcaatgctgtttcactgggtatgcggcggatccgaaaagaaaacgttgatgccggtgaacgtgcaaaaca
ggctctagcgttcgaacgactgatttcgaccagggttcgttactcatggaaaatagcgatcgctgccaggatatacgaatctggca
tttctggggattgctataacacctgttacgtatagccgaaattgccaggatcagggttaaagatatctcacgtactgacgggtggga
gaatgtaatccatattggcagaacgaaaacgctggttagcaccgagggtgtagagaaggcacttagcctgggggtaactaaactg
gtcgagcgatggatttccgtctctggtgtagctgatgatccgaaataactacctgtttgcccgggtcagaaaaatggtgttgccgcg
catctgccaccagccagctatcaactcgcgcctggaagggttttgaagcaactcatcgattgattacggcgctaaggatgactc
tggcagagatacctggcctggtctggacacagtcccgtgtcggagccgcgagatatggcccgcgtggagttcaataaccgg
agatcatgcaagctggtggctggaccaatgtaaatattgtcatgaaactatccgtaacctggatagtgaacaggggcaatggtg**c**

gctgctggaagatggcgat**taagaagttcctatactattgaaga**ataggaacttcgatcataatcagccatatcacatctgtaga
ggttttacttgctttaaaaaacctcccacacctcccctgaacctgaaacataaaatgaatgcaattggtgtgtaactggtttattgc
agcttataatggttacaataaagcaatagcatcacaattcacaataaagcattttttcactgcattctagttgtggtttgtcca
aactcatcaatgtatcttatcatgtctggatccCACTTATGCAAAGACATATATAATATATATATATATATATATATGA
TAAAGAACTTTTTTATGTTACACATTTTCCAGATATAAGAGACTGACCAGTCTTGTACAGTTTTTTTTT
TTTATTGACTGTTGGGAGTTTATGTTCTCTAGTTTTTGTGAGGTTTGACTTAATTTATATAAATACTT
TTTTTGTGTTGTTGTTTCATGTGAATGAGTGTCTAGGCAGGACCTGTGGCCAAGTTCTTAGTAGCTG
TTTATCTGTGTAGGACTGTAGAAGTGTAGAGGAAGAACTGAACATCCAGAATGTGTGGTAAAT
TGAATAAAGCTAGCCGTGATCC

Sequence highlighted in red is to be replaced by the ERT2 cassette. Nucleotides highlighted in yellow background and red are bases that are changed for silent mutation in the mutant allele.

CXCR4-CRE-ERT2-EM2-B6N

ccacgaccaagtgacagcaatgctgcttactggttatgcgcggatccgaaaagaaaacgttgatgccggtgaacgtgcaaaaca
ggctctagcgttcgaacgcactgatttcgaccagggttcactcatggaaaatagcgatcgctgccaggatatacgaatctggca
tttctggggattgcttataacaccctgttacgtatagccgaaattgccaggatcagggttaagatatctcacgtactgacggtggga
gaatgtaatccatattggcagaacgaaaacgctggttagcaccgagggtgtagagaaggcacttagcctgggggtaactaaactg
gtcagcgatggatttccgtctggtgtagctgatgatccgaataactacctgttttgcgggtcagaaaaatggtgtgccgcgc
catctgccaccagccagctatcaactcgccctggaagggattttgaagcaactcatcgattgattacggcgctaaggatgactc
tggtcagagatactggcctggtctggacacagtgcccggtcggagccgcgagataggcccgctggagtttaataccgg
agatcatgcaagctggtggctggaccaatgtaaatattgtcatgaactatatccgtaacctggatagtgaacaggggcaatggtg
AgActgctggaagatggcgat**ctcagacca**ctgctggagacatgagagctccaaccttggccaagcccgctcatgatcaaac
gctctaagaagaacagcctggccttgcctgacggcgaccagatggtcagtgcttggatgctgagccccatactctattc
cgagtatgatcctaccagacccttcagtgaagcttcgatgatgggcttactgaccaacctggcagacagggagctggttcatgat
caactgggccaagaggggtgcaggcttggattgacctccatgatcaggctcaccttctagaatgtgcttggctagagatcctg
atgattggtctcgtctggcctccatggagcaccagtgaaactactgcttcttaactgcttggacaggaaccagggaaaat
gtgtagagggcatggtggagatcttcagatgctgctggctacatcatctcggttccgcatgatgaatctgcagggagaggattgt
gtcctcaaatctatttttcttaattctggagtgtacacatttctgtccagcaccctgaagtcttggagagaaggaccatatcc
accgagtctggacaagatcacagacactttgatccacctgatggccaaggcaggcctgacctgcagcagcagcaccagcggct
ggcccagctcctctcatcctctcccacatcaggcacatgagtaacaaaggcatggagcatctgtacagcatgaagtgaagaacg
tggtgccctctatgacctgctgctggaggcggcgaccccaccgctacatgcgccactagccgtggaggggcatccgtggag
gagacggaccaaagccacttggcactgctgggcttacttcatcgattccttcaaaagtattacatcacgggggaggcagaggg
ttccctgccacagcttgatgaATAGGAACTTCgatcataatcagccatatcacatctgtagaggttttacttgctttaaaaaac
tcccacacctcccctgaaactgaaacataaaatgaatgcaattggtgtgtaactggtttattgcagcttataatggttacaataa
agcaatagcatcacaattcacaataaagcattttttcactgcattctagttgtggtttgtccaactcatcaatgtatcttatcat
gtctggatccCACTTATGCAAAGACATATATAATGATAAAGAACTTTTTTATG
TTACACATTTTCCAGATATAAGAGACTGACCAGTCTTGTACAGTTTTTTTTTTTTTTTATTGACTGTTGGGA
GTTTATGTTCTCTAGTTTTTGTGAGGTTTGACTTAATTTATATAAATACTTTTTTTTTGTTTGTGTTT

ATGTGAATGAGTGTCTAGGCAGGACCTGTGGCCAAGTTCTTAGTAGCTGTTTATCTGTGTGTAGGAC
TGTAGAACTGTAGAGGAAGAACTGAACATTCCAGAATGTGTGGTAAATTGAATAAAGCTAGCCGT
GATCC

Sequence highlighted in red and underlined is the ERT2 cassette. Sequence highlighted in red is the **stop codon**. Sequence highlighted in bold is a **linking sequence**. Nucleotides highlighted with a yellow background, in bold and red as the **silent mutations** introduced to prevent the reprocessing of the engineered allele.

Nucleotide Alignment:

```

*      20      *      40      *      60      *      80      *      100     *      120
Cxcr4_EM1 : ccacgaccaagtgcagcaatgctgtttcactggttatgcggcggatccgaaaagaaaacgttgatgccggtgaacgtgcaaacaggctctagcgttcgaacgcactgatttcgaccag : 120
Cxcr4_EM2 : ccacgaccaagtgcagcaatgctgtttcactggttatgcggcggatccgaaaagaaaacgttgatgccggtgaacgtgcaaacaggctctagcgttcgaacgcactgatttcgaccag : 120

*      140     *      160     *      180     *      200     *      220     *      240
Cxcr4_EM1 : gttcgttctactcatggaaaatagcgatcgctgccaggatatacgtaatctggcatttctggggattgcttataaacacctgttacgtatagccgaaattgccaggatcagggttaaagat : 240
Cxcr4_EM2 : gttcgttctactcatggaaaatagcgatcgctgccaggatatacgtaatctggcatttctggggattgcttataaacacctgttacgtatagccgaaattgccaggatcagggttaaagat : 240

*      260     *      280     *      300     *      320     *      340     *      360
Cxcr4_EM1 : atctcacgtactgcaggtgggagaatgttaatccatattggcagaacgaaaacgctggtagcaccgcaggtgtagagaaggcacttagcctggggtaactaaactggtcgagcgatgg : 360
Cxcr4_EM2 : atctcacgtactgcaggtgggagaatgttaatccatattggcagaacgaaaacgctggtagcaccgcaggtgtagagaaggcacttagcctggggtaactaaactggtcgagcgatgg : 360

*      380     *      400     *      420     *      440     *      460     *      480
Cxcr4_EM1 : atttcgctctggtgtagctgatgatccgaataactacctgtttgcccgggtcagaaaaaatgggttgcgcgacctgcccaccagccagctatcaactcgcgcctggaagggatt : 480
Cxcr4_EM2 : atttcgctctggtgtagctgatgatccgaataactacctgtttgcccgggtcagaaaaaatgggttgcgcgacctgcccaccagccagctatcaactcgcgcctggaagggatt : 480

*      500     *      520     *      540     *      560     *      580     *      600
Cxcr4_EM1 : tttgaagcaactcatcgattgatttacggcgctaaggatgactctggctcagagatactggcctggctctggacacagtgcccgtgctggagccgcgagatatggcccgcctggagtt : 600
Cxcr4_EM2 : tttgaagcaactcatcgattgatttacggcgctaaggatgactctggctcagagatactggcctggctctggacacagtgcccgtgctggagccgcgagatatggcccgcctggagtt : 600

*      620     *      640     *      660     *      680     *      700     *      720
Cxcr4_EM1 : tcaataccggagatcatgcaagctggtggctggaccaatgtaaatattgtcatgaactatatccgtaacctggatagtgaacaggggcaatggtgagcctgctggaagatggcgat--- : 717
Cxcr4_EM2 : tcaataccggagatcatgcaagctggtggctggaccaatgtaaatattgtcatgaactatatccgtaacctggatagtgaacaggggcaatggtgAgActgctggaagatggcgatctc : 720

*      740     *      760     *      780     *      800     *      820     *      840
Cxcr4_EM1 : ----- : -
Cxcr4_EM2 : gagccaatctgctggagacatgagagctgccaacctttggccaagcccgctcatgatcaaacgctctaaagaagaacagcctggccttgtccctgacggccgaccagatgggtcagtgccttg : 840

*      860     *      880     *      900     *      920     *      940     *      960
Cxcr4_EM1 : ----- : -
Cxcr4_EM2 : ttggatgctgagccccatactctatccgagtatgatcctaccagaccttcagtgaagcttcgatgatgggcttactgaccaacctggcagacagggagctgggtccatgatcaac : 960

*      980     *      1000    *      1020    *      1040    *      1060    *      1080
Cxcr4_EM1 : ----- : -
Cxcr4_EM2 : tggcgaagagggtgccaggcttctgtgatttgacctccatgatcaggtccacctctagaatgtgcttgcttagagatcctgatgatggctctcgtctggcctccatggagaccoca : 1080

```



```
                *      2180          *      2200          *      2220          *      2240          *      2260
Cxcr4_EM1 : CAAGTCTTAGTAGCTGTTTATCTGTGTGTAGGACTGTAGAACTGTAGAGGAAGAAACTGAACATCCAGAATGTGTGGTAAATTGAATAAAGCTAGCCGTGATCC : 1335
Cxcr4_EM2 : CAAGTCTTAGTAGCTGTTTATCTGTGTGTAGGACTGTAGAACTGTAGAGGAAGAAACTGAACATCCAGAATGTGTGGTAAATTGAATAAAGCTAGCCGTGATCC : 2266
```


QC strategy employed at Harwell to check the edited allele:

Genomic DNA was extracted from ear clip biopsies and amplified in a PCR reaction using the following conditions/primer sequences:

Geno_Cxcr4_CRE-ERT2_F6 (5'-3')	CCACGACCAAGTGACAGCAAT
Geno_Cxcr4_CRE-ERT2_R6 (5'-3')	GGATCACGGCTAGCTTTATTCA
Taq Polymerase used	Roche KAPA Hotstart PCR Kit
Annealing Temperature (°C)	56
Elongation time (min)	2
WT product size (bp)	1335
Mutant product size (bp)	2266
Notes	This PCR sits externally to the donor template applied.

Geno_Cxcr4_CRE-ERT2_F3 (5'-3')	TCAAACGCTCTAAGAAGAACAGC
Geno_Cxcr4_CRE-ERT2_R7 (5'-3')	GATGTGGGAGAGGATGAGGA
Taq Polymerase used	ThermoFisher SuperFi II PCR Kit
Annealing Temperature (°C)	60
Elongation time (min)	0.5
WT product size (bp)	N/A
Mutant product size (bp)	650
Notes	This PCR is used to screen for the presence of the ERT2 cassette.

Cre_F4 (5'-3')	tccgtctctggtgtagctga
Geno_Cxcr4_CRE-ERT2_R7 (5'-3')	GATGTGGGAGAGGATGAGGA
Taq Polymerase used	ThermoFisher SuperFi II PCR Kit
Annealing Temperature (°C)	60
Elongation time (min)	0.75
WT product size (bp)	N/A

Mutant product size (bp)	1063
Notes	This PCR is used to screen whether the 5' end of the donor has gone in on-target.

Geno_Cxcr4_CRE-ERT2_F8 (5'-3')	ggctacatcatctcggttcc
Geno_Cxcr4_CRE-ERT2_R8 (5'-3')	ATAAGTGggatccagacatg
Taq Polymerase used	ThermoFisher SuperFi II PCR Kit
Annealing Temperature (°C)	60
Elongation time (min)	1
WT product size (bp)	N/A
Mutant product size (bp)	765
Notes	This PCR is used to screen whether the 3' end of the donor has gone in on-target.

All amplicons were sent for Sanger sequencing to check for integration of the donor oligo sequence at the target site. F1 sequences should be heterozygous unless on sex chromosome.

Additional integrations of the donor sequence

Copy counting of the donor sequence was carried out by ddPCR at the F1 stage to confirm donor oligos were inserted once on target into the genome. The following Taqman assay was used to copy count the donor sequence compared against a VIC-labelled reference assay for Dot1l:

Assay name	ERT2-MUT1
Forward Primer (5'-3')	GCTTCGATGATGGGCTTACT
Reverse Primer (5'-3')	CCTGATCATGGAGGGTCAAATC
Probe (5'-3')	TGGTTCACATGATCAACTGGGCGA
Label	FAM

This ddPCR assay is specific to the donor used to create the engineered mutation and only mutant alleles are expected to be recognised by this assay. Therefore, WT controls are expected to call at 0 copies and a single integration for a correct mutation is expected to call at 1 copy for F1 (HET) animals.

Reference Assay Name	Dot1l
Forward primer (5'-3')	GCCCCAGCACGACCATT
Reverse primer (5'-3')	TAGTTGGCATCCTTATGCTTCATC
Probe (5'-3')	CCCAACAGGCCTGGATTCTCAATGC
Label	VIC

VIC-labelled reference assay for Dot1l gene.

No additional donor integrations were detected in the animals taken forward to establish the colony.