

PCR TTR cre tam
-----------------

## 1- PCR mix

2,50 µl 10X buffer(Q BIOgene ; 1,5 mM MgCl<sub>2</sub>)  
 0,25 µl dNTP (20 mM total)  
 0,15 µl primer cre1 (200 ng/µl)  
 0,15 µl primer cre22 (200 ng/µl)  
 0,20 µl Taq polymerase (Q BIOgene ; 5U/µl)  
 19,75 µl distilled water

## 2- PCR reaction :

23,0 µl mix  
 2,0 µl DNA (200ng/µl)

## 3- PCR program :

5 min	94°C	
1 min	94°C	
1 min	57°C	35 cycles
2 min	72°C	
5 min	72°C	

## 4- Primers:

- Primer cre1                      5' CTG CAT TAC CGG TCG ATG CA 3'  
 - Primer cre2                      5' ACG TTC ACC GGC ATC AAC GT 3'

## 5- Amplifications:

transgene :                      380 bp