

STRAIN NAME:**NGF 13 KO****REAGENTS AND SOLUTIONS:**

- Taq polymerase 10X Buffer (1,5 mM MgCl₂) (MP QBiogene)
- dNTP (MPQ Biogene)
- Primers (Invitrogen)
- Taq polymerase (MPQBiogene)
- Mineral oil (EUROBIO)

REACTION MIX:

- Taq Buffer 1X
- dNTP (0,4 mM total)
- primer NGF13F (0,36 µM)
- primer NGF13R (0,72 µM)
- primer NGF13KO (0,36 µM)
- Taq polymerase (40U/ml)
- DNA

PROGRAM:

5 min	94°C	35 cycles
30 sec	94°C	
30 sec	60°C	
1 min	72°C	
5 min	72°C	

MIGRATION :

Migration of the samples on a 1.5% agarose gel.

PRIMERS SEQUENCES

- | | |
|------------------|-------------------------------------|
| - Primer NGF13F | 5' TCTCTGAACAGATTTCTGGGAGAGGTAGC 3' |
| - Primer NGF13R | 5' CAAGATTTACCACCAACCCAACTGTGC 3' |
| - Primer NGF13KO | 5' GCTAGGGAAGGAGAAGAATGGGAACG3' |

AMPLIFICATION

Mutant allele : 700 bp
WT allele :415 bp