



Genotyping protocol

Cdk8

IR00003723 / G13

(ICS internal reference)

This report has been prepared by: **Pauline Cayrou**
33 (0)3 88 65 56 55
genotyping@igbmc.fr

This report has been validated by: **Sylvie Jacquot, PhD, Head of Genotyping Service**
33 (0)3 88 65 57 44
genotyping@igbmc.fr

The first version of this report was finalized the: 14 Oct 2013

For any question, please contact:

Institut Clinique de la Souris - ICS - Mouse Clinical Institute
1 rue Laurent Fries, BP 10142
67404 Illkirch Cedex, France
Email: genotyping@igbmc.fr
Web site: <http://www-mci.u-strasbg.fr/>

TABLE OF CONTENTS

Table of contents2

1. Genotyping protocol and data2

 1.1. Genotyping strategy2

 1.2. PCR protocol4

 1.3. Picture of genotyping with various alleles5

2. Cre and Flp genotyping method6

1. Genotyping protocol and data


This section describes the condition used at the Mouse Clinical Institute (ICS) to genotype your **Cdk8** Constitutive Knockout / Conditional Knockout (KO-cKO pc) project.

1.1. Genotyping strategy

The map below describes the position of the primers used for genotyping for each possible allele.



KO-cKO pc _ Genotyping strategy



Sequence of primers used for genotyping:

Position	Primers	Sequence
Lxr	3255	ACTGATGGCGAGCTCAGACCATAAC
Kr	3278	GGGCAAGAACATAAAGTGACCCTCC
Rif	5966	GCACATGGCTGAATATCGACGGT
Ef	6108	GCTTCCAAAAAAGAGTACCAAAGGCG
Ef ²	6109	CGTAGGTAGCAATCTGGTCGGGGT
L3f	6110	TTATGGGTGAAGGTCCTGAGCGA
L3f ²	6111	ACCTCCCCCTGTCTCAGACAGTCCTT
L3r	6112	GGAAAACGACAATGGAAGCAGCAG
Er	6113	CAGGTACACAGGCTGGATTTGCAC

². For a selected position, a second primer was designed

PCR fragments expected size (bp):

Region analyzed	Primers used	Position on the primer (see the map above)	Targeted allele (HR)	cKO allele	KI allele	WildType allele (WT)
5' part of the selection marker	6108-3278	Ef / Kr	262	---	---	---
Presence of the distal loxP	6111-6112	L3f ² / L3r	499	499	---	446
Distal loxP specific PCR	6110-3255	L3f / Lxr	472	472	---	---
Excision of the selection marker	6109-6113	Ef ² / Er	7501*	597	---	414
Cre total excision	5966-3255	Ri1f / Lxr	3472*	---	471**	---

*: this PCR product will not be observed using our PCR genotyping conditions (see description below)

** : this PCR is only verified if mice are generated

---: no Amplicon should be obtained

1.2. PCR protocol

This section describes the composition of the mix and cycling conditions used for genotyping.

Reagents:	Volume:
- FastStart PCR Master (Roche)	7.5µl
- DNA (50ng/µl)	1.5µl
- 5' primer (100 µM)	0.06µl
- 3' primer (100 µM)	0.06µl
- Sterile H ₂ O	up to 15 µl

Cycling conditions:


Temp	Time	#Cycles
95°C	4min	1
94°C	30s	34
62°C	30s	
72°C	1min	
72°C	7min	1
20°C	5min	1

NB: These PCR conditions have been optimized for high-throughput genotyping. Adaptation to small-scale may be required.

1.3. Picture of genotyping with various alleles

Analysis of PCR products pattern was not done by gel electrophoresis but using LabChip® 90 microfluidic apparatus. PCR products were run on the HT DNA 5K LabChip® 90 Assay Kit.

Representative genotyping picture



Note that as this technology is more sensitive than gel analysis, non specific signals and/or primer dimers may be visible on the picture.

2. Cre and Flp genotyping method

You will find the genotyping protocol in the publication:

[Highly-efficient, fluorescent, locus directed cre and FlpO deleter mice on a pure C57BL/6N genetic background.](#)

Birling MC, Dierich A, Jacquot S, Hérault Y, Pavlovic G.

Genesis. 2012 Jun;50(6):482-9. doi: 10.1002/dvg.20826. Epub 2012 Mar 20.