



STRAIN NAME :

EGFP

REAGENTS AND SOLUTIONS:

- Kapa Mouse Genotyping kit (Sigma)
- Primers 30µM (Eurofins)
- Mineral Oil (Fisher Scientific)
- DNA (extracted with Kapa Mouse Genotyping kit)

REACTION MIX :

MIX	Quantity / well
Water	6,80 µl
2X Kapa Buffer	10,00 µl
EGF1	0,60 µl
EGF2	0,60 µl
DNA	2,00 µl
Final Quantity	20,00 µl

PCR PROGRAM:

3 min	95°C	
15 sec	94°C	
15 sec	60°C	35 cycles
30 sec	72°C	
5 min	72°C	

MIGRATION :

Migration on 1.5% agarose gel.

PRIMERS SEQUENCES

- EGF1 5' TGAACCGCATCGAGCTGAAGGG 3'
- EGF2 5' TCCAGCAGGACCATGTGATCGC 3'

EXPECTED AMPLIFICATIONS

EGFP : 307 bp

Perfection date : 25/11/2013

Update :