

STRAIN NAME :

**HSD17B4 Flox**

REAGENTS AND SOLUTIONS:

- Kapa Mouse Genotyping kit (Sigma)
- Primers 30µM (Eurofins)
- Mineral Oil (Fisher Scientific)
- DNA (extracted with Kapa Mouse Genotyping kit)

REACTION MIX :

MIX	Quantity / well	
Water	6,80	µl
2X Kapa Buffer	10,00	µl
HSD17B4FloxF	0,60	µl
HSD17B4FloxR	0,60	µl
DNA	2,00	µl
Final Quantity	20,00	µl

PCR PROGRAM:

3 min	95°C	
15 sec	94°C	
15 sec	60°C	35 cycles
30 sec	72°C	
5 min	72°C	

MIGRATION :

Migration on 1.5% agarose gel.

PRIMERS SEQUENCES

- HSD17B4FloxF      5' CCCAACGCTGGGTCACGGATGACG 3'
- HSD17B4FloxR    5' GCAACCATAAGTTACACAAAATGCC 3'

EXPECTED AMPLIFICATIONS

WT allele : 256 bp  
 Modified allele: 304 bp

*Perfection date : 03/12/2021*

*Update :*