

STRAIN NAME :

Wls cKO

REAGENTS AND SOLUTIONS:

- Kapa Mouse Genotyping kit (Sigma)
- Primers 30µM (Eurofins)
- Mineral Oil (Fisher Scientific)
- DNA (extracted with Kapa Mouse Genotyping kit)

REACTION MIX :

MIX	Quantity / well
Water	6,80 µl
2X Kapa Buffer	10,00 µl
WlsF	0,60 µl
WlsR	0,60 µl
DNA	2,00 µl
Final Quantity	20,00 µl

PCR PROGRAM:

3 min	95°C	
15 sec	94°C	35 cycles
15 sec	60°C	
30 sec	72°C	
5 min	72°C	

MIGRATION :

Migration on 1.5% agarose gel.

PRIMERS SEQUENCES

- WlsF 5' CCCCTTTCCCTCTCGGTTCC 3'
- WlsR 5' GCGGCATGGAAGCCAAGGGC 3'

EXPECTED AMPLIFICATIONS

WT allele (+) : 267 bp
Flox allele (f) : 353 bp

Validation date : 21/01/2026

Update :