

PCR CD11cDTR

1- PCR mix

2,50 µl 10X buffer(Q BIOgene ; 1,5 mM MgCl₂)
 0,25 µl dNTP (20 mM total)
 0,20 µl primer DTR1 (200 ng/µl)
 0,20 µl primer DTR2 (200 ng/µl)
 0,20 µl Taq polymerase (Q BIOgene ; 5U/µl)
 19,65 µl distilled water

2- PCR reaction :

23,0 µl mix
 2,0 µl DNA (200ng/µl)

3- PCR program :

5 min	94°C	
1 min	94°C	
1 min	60 °C	35 cycles
2 min	72°C	
5 min	72°C	

4- Primers:

- Primer DTR1 5' GCC ACC ATG AAG CTG CTG CCG 3'
 - Primer DTR2 5' TCA GTG GGA ATT AGT CAT GCC 3'

5- Amplification:

transgene : 625 bp