

STRAIN: EM04733

PRIMER mapk8ip for 5'-CGG CAG AGA GTG AAA TCA AGC TTG-3'

PRIMER mapk8ip rev 5'-CTC CTA CAC AGA AAC AGG CCT ACA-3'

PRIMER FRT 24 REV 5'-CCT TCC TCC TAC ATA GTT GGC AGT-3'

PRIMER mapk8ip for 15pm / total volume 30ul Tm: 72°C

PRIMER mapk8ip rev 15 pm / total volume 30ul Tm: 72°C

PRIMER FRT 24 REV 15pm / total volume 30ul Tm: 72°C

PCR MIX

DNTPS (2mM) 3ul

PCR BUFFER 10x 3ul

PRIMER mapk8ip for (0.5pm/ul)

PRIMER mapk8ip rev (0.5pm/ul)

PRIMER FRT 24 REV (0.5pm/ul)

H2O

TAQ POLYMERASE (5U/ul) 0.3ul

DNA 2ul

TOTAL 30ul

PCR PROGRAM:

- 1) 95°C x 2'00''
- 2) 95°C x 0'30''
- 3) 46°C x 0'30''
- 4) 72°C x 1'00''
- 5) step 2 to step 5 x 50 repeats
- 6) 72°C x 10'00''
- 7) 4°C forever
- 8) End

MUTATION: 190 BP

WILD TYPE: 110 BP